Mibefradil

Cat. No.:	HY-15553	
CAS No.:	116644-53-2	
Molecular Formula:	C ₂₉ H ₃₈ FN ₃ O ₃	
Molecular Weight:	495.63	
Target:	Calcium Channel	
Pathway:	Membrane Transporter/Ion Channel; Neuronal Signaling	
Storage:	Please store the product under the recommended conditions in the Certificate of Analysis.	

SIDEOGICAE ACTIVITY		
Description	Mibefradil (Ro 40-5967) is a calcium channel blocker with moderate selectivity for T-type Ca ²⁺ channels displaying IC ₅₀ s of 2.7 μM and 18.6 μM for T-type and L-type currents, respectively ^[1] .	
IC ₅₀ & Target	IC50: 2.7 μM (T-type calcium channel), 18.6 μM (L-type calcium channel) $^{[1]}$	
In Vitro	Mibefradil inhibits reversibly the T- and L-type currents with IC ₅₀ values of 2.7 and 18.6 μM, respectively. The inhibition of the L-type current is voltage-dependent, whereas that of the T-type current is not. Ro 40-5967 blocks T-type current already at a holding potential of -100 mV ^[1] At a higher concentration (20 μM), Mibefradil reduces the amplitude of excitatory junction potentials (by 37±10 %), slows the rate of repolarisation (by 44±16 %) and causes a significant membrane potential depolarisation (from -83±1 mV to -71±5 mV). At a higher Mibefradil concentration (20 μM) there is significant membrane potential depolarisation and a slowing of repolarisation. These actions of Mibefradil are consistent with K ⁺ channel inhibition, which has been shown to occur in human myoblasts and other cells ^[2] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.	
In Vivo	The hearing thresholds of the 24-26 week old C57BL/6J mice differed following the 4-week treatment period. The hearing threshold at 24 kHz is significantly decreased in the Mibefradil-treated and benidipine-treated groups compared with the saline-treated group, rats receiving Mibefradil or Ethosuximide show	

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

significant lower Ca_V3.2 expression in the spinal cord and DRG^[4].

PROTOCOL

Animal Administration ^{[3][4]}

Mice^[3]

A total of 30 male C57BL/6J mice (age, 6-8 weeks) are randomized into three groups for the detection of three calcium channel receptor subunits α1G, α1H and α1I, using reverse transcription-quantitative polymerase chain reaction (RT-qPCR). In addition, a further 30 C57BL/6J male mice (age, 24-26 weeks) are allocated at random into three treatment groups: Saline, Mibefradil and benidipine. Each group is subjected to auditory brainstem recording (ABR) and distortion product otoacoustic emission (DPOAE) tests following treatment. Mibefradil and benidipine are dissolved in physiological saline solution. A preliminary experiment led to the selection of dosages of 30 mg/kg/day Mibefradil and 10 mg/kg/day Benidipine. The drugs are administered to the mice by gavage for four consecutive weeks. Rats^[4]

Male Sprague-Dawley rats (200-250 g) are used for right L5/6 SNL to induce neuropathic pain. Intrathecal infusion of saline or TCC blockers [Mibefradil (0.7 μg/h) or Ethosuximide (60 μg/h)] is started after surgery for 7 days. Fluorescent immunohistochemistry and Western blotting are used to determine the expression pattern and protein level of Ca_V3.2. Hematoxylin-eosin and toluidine blue staining are used to evaluate the neurotoxicity of tested agents. MCE has not independently confirmed the accuracy of these methods. They are for reference only.

CUSTOMER VALIDATION

- Br J Pharmacol. 2021 Jan;178(2):346-362.
- J Cell Physiol. 2021 Mar 11.
- Front Pharmacol. 23 February 2022.
- Mediat Inflamm. 2020 Nov 10;2020:3691701.
- Eur J Pharmacol. 2021 Feb 5;892:173782.

See more customer validations on www.MedChemExpress.com

REFERENCES

[1]. Mehrke G, et al. The Ca(++)-channel blocker Ro 40-5967 blocks differently T-type and L-type Ca++ channels. J Pharmacol Exp Ther. 1994 Dec;271(3):1483-8.

[2]. Brain KL, et al. The sources and sequestration of Ca(2+) contributing to neuroeffector Ca(2+) transients in the mouse vas deferens. J Physiol. 2003 Dec 1;553(Pt 2):627-35.

[3]. Yu YF, et al. Protection of the cochlear hair cells in adult C57BL/6J mice by T-type calcium channel blockers. Exp Ther Med. 2016 Mar;11(3):1039-1044.

[4]. Shiue SJ, et al. Chronic intrathecal infusion of T-type calcium channel blockers attenuates CaV3.2 upregulation in nerve-ligated rats. Acta Anaesthesiol Taiwan. 2016 Oct 17. pii: S1875-4597(16)30071-6.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898 Fax: 609-228-5909 E-mail: tech@MedChemExpress.com Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA