

Basement Membrane Matrix HC

1 Contents

Components	HY-K6008-5 mL	HY-K6008-10 mL
Basement Membrane Matrix HC	5 mL	10 mL

2 Introduction

MedChemExpress Basement Membrane Matrix HC is a natural basement membrane matrix extracted from mouse tumors and is composed mainly of various growth factors and extracellular matrix components. The main extracellular matrix components are: Laminin, Col-IV, Entactin, Heparan sulphate proteoglycans, etc. The main growth factor components are: insulin-like growth factor (IGF-1), transforming growth factor beta (TGF- β), vascular endothelial growth factor (VEGF), epidermal growth factor (EGF), fibroblast growth factor (bFGF), etc. The percentage of each of these growth factors is higher than that of standard basement membrane matrix.

MedChemExpress Basement Membrane Matrix HC is mainly used for angiogenesis experiments, animal models and 3D tumor model construction.

3 Characteristics

Source	Mouse Tumor
Color	Color Yellow
Appearance	Appearance $\leq 0^{\circ}\text{C}$:Solid; $0-4^{\circ}\text{C}$:Liquid; $\geq 4^{\circ}\text{C}$:Semi-gel or Gel
Protein concentration	Protein concentration 16~26 mg/mL
Endotoxin	Endotoxin ≤ 4.5 EU/ mL
Gelling time	Gelling time Room temperature: 5-30 min
Formulation	Formulation Supplied in DMEM and 50 $\mu\text{g}/\text{mL}$ gentamicin

4 Storage

Store at -20°C , 2 years.

Avoid repeated freezing and thawing.

5 Precautions

1. Thaw matrix with packaging in a 4°C refrigerator. After thawing, make aliquots and keep them frozen.
2. Use pre-cool consumables to avoid gelling.
3. Avoid holding the container in your hands which may cause semi-gel. If it happens, please put matrix back to 0-4°C refrigerator for 1-2 hours to restore its fluidity.
4. Color variation is normal and does not affect the use of matrix. Due to the interaction of carbon dioxide with bicarbonate buffer and phenol red, the color may vary from straw yellow to deep red.
5. This product is for R&D use only, not for drug, household, or other uses.
6. For your safety and health, please wear a lab coat and disposable gloves to operate.