


Omipalisib

Cat. No.:	HY-10297		
CAS No.:	1086062-66-9		
Molecular Formula:	C ₂₅ H ₁₇ F ₂ N ₅ O ₃ S		
Molecular Weight:	505.5		
Target:	PI3K; mTOR; Autophagy		
Pathway:	PI3K/Akt/mTOR; Autophagy		
Storage:	Powder	-20°C	3 years
		4°C	2 years
	In solvent	-80°C	6 months
		-20°C	1 month



SOLVENT & SOLUBILITY

In Vitro	DMSO : 50 mg/mL (98.91 mM; Need ultrasonic)				
		Solvent Concentration	Mass		
	Preparing Stock Solutions		1 mg	5 mg	10 mg
		1 mM	1.9782 mL	9.8912 mL	19.7824 mL
5 mM		0.3956 mL	1.9782 mL	3.9565 mL	
	10 mM	0.1978 mL	0.9891 mL	1.9782 mL	
Please refer to the solubility information to select the appropriate solvent.					
In Vivo	1. Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline Solubility: ≥ 2.5 mg/mL (4.95 mM); Clear solution				

BIOLOGICAL ACTIVITY

Description	Omipalisib (GSK2126458) is an orally active and highly selective inhibitor of PI3K with K _i s of 0.019 nM/0.13 nM/0.024 nM/0.06 nM and 0.18 nM/0.3 nM for p110α/β/δ/γ, mTORC1/2, respectively. Omipalisib has anti-cancer activity ^{[1][2][3]} .			
IC₅₀ & Target	p110α 0.019 nM (K _i)	p110α-E545K 0.008 nM (K _i)	p110α-E542K 0.008 nM (K _i)	p110α-H1047R 0.009 nM (K _i)
	p110β 0.13 nM (K _i)	p110δ 0.024 nM (K _i)	p110γ 0.06 nM (K _i)	mTORC1 0.18 nM (K _i)
	mTORC2 0.3 nM (K _i)			

In Vitro	<p>Omipalisib (GSK2126458) potently inhibits the activity of common activating mutants of p110α (E542K, E545K, and H1047R) found in human cancer with K_i of 8 pM, 8 pM and 9 pM, respectively. Omipalisib causes a significant reduction in the levels of pAkt-S473 with remarkable potency in T47D and BT474 cells with IC₅₀ of 0.41 nM and 0.18 nM, respectively. Furthermore, Omipalisib (GSK2126458) leads to a G1 cell cycle arrest and produces the inhibitory effect on cell proliferation in a large panel of cell lines, including T47D and BT474 breast cancer lines with IC₅₀ of 3 nM and 2.4 nM, respectively^[1]. The combination of Omipalisib or GSK1120212 with Omipalisib enhances cell growth inhibition and decreases S6 ribosomal protein phosphorylation in drug-resistant clones from the A375 BRAF(V600E) and the YUSIT1 BRAF(V600K) melanoma cell lines^[2]. Omipalisib (GSK2126458) potentiates the antiproliferative activity of DDR1-IN-1 in colorectal cancer cell lines^[3]. MCE has not independently confirmed the accuracy of these methods. They are for reference only.</p>
In Vivo	<p>In a BT474 human tumor xenograft model, Omipalisib (GSK2126458) treatment results in a dose-dependent reduction in pAkt-S473 levels, and exhibits dose-dependent tumor growth inhibition at a low dose of 300 μg/kg. Besides, Omipalisib (GSK2126458) shows low blood clearance and good oral bioavailability in four preclinical species (mouse, rat, dog, and monkey)^[1]. MCE has not independently confirmed the accuracy of these methods. They are for reference only.</p>

PROTOCOL

Cell Assay ^[1]

BT474, HCC1954 and T-47D (human breast) are cultured in RPMI-1640 containing 10% fetal bovine serum at 37°C in 5% CO₂ incubator. Cells are split into T75 flask two to three days prior to assay set up at density which yields approximately 70-80% confluence at time of harvest for assay. Cells are harvested using 0.25% trypsin-EDTA. Cell counts are performed on cell suspension using Trypan Blue exclusion staining. Cells are then plated in 384 well black flat bottom polystyrene in 48 μ L of culture media per well at 1,000 cells/well. All plates are placed at 5% CO₂, 37°C overnight and Omipalisib (GSK2126458) is added the following day. One plate is treated with CellTiter-Glo for a day 0 (t=0) measurement and read as described below. Omipalisib (GSK2126458) is prepared in clear bottom polypropylene 384 well plates with consecutive two fold dilutions. 4 μ L of these dilutions are added to 105 μ L culture media, after mixing the solution, 2 μ L of these dilutions are added into each well of the cell plates. The final concentration of DMSO in all wells is 0.15%. Cells are incubated at 37°C, 5% CO₂ for 72 hours. Following 72 hours of incubation with Omipalisib each plate is developed and read. CellTiter-Glo reagent is added to assay plates using a volume equivalent to the cell culture volume in the wells. Plates are shaken for approximately two minutes and incubated at room temperature for approximately 30 minutes and chemiluminescent signal is read on the Analyst GT reader. Results are expressed as a percent of the t=0 and plotted against the Omipalisib (GSK2126458) concentration. Cell growth inhibition is determined for Omipalisib (GSK2126458) by fitting the dose response with a 4 or 6 parameter curve fit using XLfit software and determining the concentration that inhibits 50% of the cell growth (GI₅₀) with the Y min as the t=0 and Y max as the DMSO control. Value from wells with no cells is subtracted from all samples for background correction.. MCE has not independently confirmed the accuracy of these methods. They are for reference only.

CUSTOMER VALIDATION

- J Hepatol. 2021 Aug;75(2):363-376.
- Sci Transl Med. 2018 Jul 18;10(450):eaq1093.
- Clin Cancer Res. 2014 Nov 1;20(21):5483-95.
- Cell Syst. 2020 Jan 22;10(1):66-81.e11.
- Cell Syst. 2020 Jan 22;10(1):66-81.e11.

See more customer validations on www.MedChemExpress.com

REFERENCES

-
- [1]. Knight SD, et al. Discovery of GSK2126458, a Highly Potent Inhibitor of PI3K and the Mammalian Target of Rapamycin. ACS Med. Chem. Lett. 2010, 1 (1), 39-43.
- [2]. Greger JG, et al. Combinations of BRAF, MEK, and PI3K/mTOR inhibitors overcome acquired resistance to the BRAF inhibitor GSK2118436 dabrafenib, mediated by NRAS or MEK mutations. Mol Cancer Ther. 2012 Apr;11(4):909-20.
- [3]. Kim HG, et al. Discovery of a potent and selective DDR1 receptor tyrosine kinase inhibitor. ACS Chem Biol. 2013 Oct 18;8(10):2145-50.
-

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA