

DEAE Cross-linked dextran A 25

Cat. No.:	HY-103703
CAS No.:	12609-80-2
Target:	Others
Pathway:	Others
Storage:	4°C, protect from light * In solvent : -80°C, 6 months; -20°C, 1 month (protect from light)

DEAE Cross-linked dextran A 25

SOLVENT & SOLUBILITY

In Vitro	DMF : < 1 mg/mL (ultrasonic;warming;heat to 60°C) (insoluble) DMSO : < 1 mg/mL (ultrasonic;warming;adjust pH to 5 with HCl;heat to 60°C) (insoluble or slightly soluble)
----------	---

BIOLOGICAL ACTIVITY

Description	DEAE Cross-linked dextran A 25 is a weak anion exchanger ^[1] .
-------------	---

REFERENCES

[1]. Masataka Hiraide, et al. Evaluation of humic complexes of trace metals in river water by adsorption on indium-treated XAD-2 resin and DEAE-Sephadex A-25 anion exchanger. Fresenius' Journal of Analytical Chemistry volume 348, pages758–761 (1994).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA