


BODIPY TMR-X SE

Cat. No.:	HY-110211
CAS No.:	217190-15-3
Molecular Formula:	C ₃₁ H ₃₅ BF ₂ N ₄ O ₆
Molecular Weight:	608.44
Target:	Fluorescent Dye
Pathway:	Others
Storage:	Please store the product under the recommended conditions in the Certificate of Analysis.



BIOLOGICAL ACTIVITY

Description

BODIPY TMR-X SE is a potent fluorescent dye. BODIPY TMR-X SE can be used to label the primary amines (R-NH₂) of proteins, amine-modified oligonucleotides, and other amine-containing molecules. BODIPY TR-X binds to protein or antibody and has bright, orange fluorescent light. (λ_{ex} =544 nm, λ_{em} =570 nm)^{[1][2]}.

REFERENCES

[1]. Boersma MD, et, al. Hydrophile scanning as a complement to alanine scanning for exploring and manipulating protein-protein recognition: application to the Bim BH3 domain. Protein Sci. 2008 Jul;17(7):1232-40.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA