


IPSU

Cat. No.:	HY-13796		
CAS No.:	1373765-19-5		
Molecular Formula:	C ₂₃ H ₂₇ N ₅ O ₂		
Molecular Weight:	405.49		
Target:	Orexin Receptor (OX Receptor)		
Pathway:	GPCR/G Protein; Neuronal Signaling		
Storage:	Powder	-20°C	3 years
		4°C	2 years
	In solvent	-80°C	2 years
		-20°C	1 year



SOLVENT & SOLUBILITY

In Vitro

DMSO : ≥ 30 mg/mL (73.98 mM)
 * "≥" means soluble, but saturation unknown.

Preparing Stock Solutions	Solvent	Mass	1 mg	5 mg	10 mg
	Concentration				
	1 mM		2.4662 mL	12.3308 mL	24.6615 mL
	5 mM		0.4932 mL	2.4662 mL	4.9323 mL
	10 mM		0.2466 mL	1.2331 mL	2.4662 mL

Please refer to the solubility information to select the appropriate solvent.

BIOLOGICAL ACTIVITY

Description	IPSU is a selective, orally available and brain penetrant OX2R antagonist with a pK _i of 7.85.
IC₅₀ & Target	pK _i : 7.85 (OX2R), 6.29 (OX1R) ^[1]
In Vitro	Orexin receptor antagonists represent attractive targets for the development of drugs for the treatment of insomnia. IPSU binds rapidly and reaches equilibrium very quickly in binding and/or functional assays ^[2] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.
In Vivo	IPSU has low blood clearance, shows high maximal blood exposure and AUC after oral dosing. It exhibits an acceptable absolute oral bioavailability and a brain/blood concentration ratio that indicated favorable brain penetration. IPSU increases sleep when dosed during the mouse active phase (lights off); IPSU induces sleep primarily by increasing NREM sleep. IPSU shows a fast onset of action, with a clear increase in total sleep time during the first hour after dosing. The effect lasts 4-5 h, after which time the total sleep time per hour is the same as on vehicle day ^[1] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.

PROTOCOL

Kinase Assay ^[2]

Competition experiments are performed with a single concentration of radioligand and six concentrations of competitor (unlabeled ligands; BBAC, almorexant, SB-649868, suvorexant, filorexant or IPSU). 4.6 nM [³H]-BBAC is added simultaneously with various concentrations of unlabeled ligand (0.1 nM-10 μM) to membranes (150 μL/well) in 50 μL/well of assay buffer with a total volume of 250 μL/well. The amount of [³H]-BBAC bound to receptors is determined at room temperature at different time points (ranging from 15 min to 4 h) and terminated by rapid vacuum filtration and liquid scintillation counting^[2].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

Animal Administration ^[1]

Mice: Freely moving C57Bl/6 mice with chronically implanted electrodes are well habituated to the experiment boxes and had access to food and water ad libitum. The test compounds (IPSU) or vehicle are administered per os as a suspension in 0.5% methylcellulose immediately prior to lights off and start of recording. Movement is recorded using infrared sensors in the roof of the box. EEG/EMG signals and motility data are used to score 10 s epochs into wake, NREM sleep, and REM sleep. Each animal served as its own control by application and recording of vehicle the day before compound (IPSU) dosing^[1]. MCE has not independently confirmed the accuracy of these methods. They are for reference only.

CUSTOMER VALIDATION

- Molecules. 2020 Feb 25;25(5):1018.

See more customer validations on www.MedChemExpress.com

REFERENCES

[1]. Betschart C, et al. Identification of a novel series of orexin receptor antagonists with a distinct effect on sleep architecture for the treatment of insomnia. J Med Chem. 2013 Oct 10;56(19):7590-607.

[2]. Callander GE, et al. Kinetic properties of "dual" orexin receptor antagonists at OX1R and OX2R orexin receptors. Front Neurosci. 2013 Dec 3;7:230.

[3]. Hoyer D, et al. Distinct effects of IPSU and suvorexant on mouse sleep architecture.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA