CamA-IN-1

Cat. No.:	HY-148578			
CAS No.:	2924055-46-7			
Molecular Formula:	$C_{23}H_{36}N_6O_6$			
Molecular Weight:	492.57			
Target:	Bacterial			
Pathway:	Anti-infection			
Storage:	Powder	-20°C	3 years	
		4°C	2 years	
	In solvent	-80°C	6 months	
		-20°C	1 month	

®

MedChemExpress

SOLVENT & SOLUBILITY

In Vitro	DMSO : 100 mg/mL (203.02 mM; Need ultrasonic)					
	Solvent Mass Concentration	1 mg	5 mg	10 mg		
	Preparing Stock Solutions	1 mM	2.0302 mL	10.1508 mL	20.3017 mL	
		5 mM	0.4060 mL	2.0302 mL	4.0603 mL	
	10 mM	0.2030 mL	1.0151 mL	2.0302 mL		
	Please refer to the so	lubility information to select the app	propriate solvent.			
In Vivo		one by one: 10% DMSO >> 40% PEC (mL (5.08 mM); Clear solution; Need		0 >> 45% saline		
2. Add each solvent one by one: 10% DMSO >> 90% (20% SBE-β-CD in saline) Solubility: 2.5 mg/mL (5.08 mM); Clear solution; Need ultrasonic						
		one by one: 10% DMSO >> 90% cor ′mL (5.08 mM); Clear solution; Need				

BIOLOGICAL ACTIV	
Description	CamA-IN-1 is a <i>Clostridioides difficile</i> -specific DNA adenine methyltransferase (CamA) inhibitor. CamA-IN-1 has inhibitory for CamA with IC ₅₀ and K _d values of 0.4 μM and 0.2 μM, respectively. CamA-IN-1 can be used for the research of Clostridioides difficile infections (CDI) ^[1] .
IC ₅₀ & Target	IC50: 0.4 μM (CamA) ^[1] Kd: 0.2 μM (CamA) ^[1]

H

ΗŅ

In Vitro

CamA-IN-1 (compound 39) has inhibitory for CamA with IC₅₀ and K_d values of 0.4 µM and 0.2 µM, respectively^[1]. CamA-IN-1 is selective for CamA over closely related bacterial and mammalian DNA and RNA MTases, protein lysine and arginine MTases, and human adenosine receptors^[1].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

REFERENCES

[1]. Jujun Zhou, et al. Systematic Design of Adenosine Analogs as Inhibitors of a Clostridioides difficile- Specific DNA Adenine Methyltransferase Required for Normal Sporulation and Persistence. J Med Chem. 2023 Jan 12;66(1):934-950.

Caution: Product has not been fully validated for medical applications. For research use only.

 Tel: 609-228-6898
 Fax: 609-228-5909
 E-mail: tech@MedChemExpress.com

 Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA