


GKA50

Cat. No.:	HY-15671	
CAS No.:	851884-87-2	
Molecular Formula:	C ₂₆ H ₂₈ N ₂ O ₆	
Molecular Weight:	464.51	
Target:	Glucokinase	
Pathway:	Metabolic Enzyme/Protease	
Storage:	Powder	-20°C 3 years
	In solvent	-80°C 6 months
		-20°C 1 month



SOLVENT & SOLUBILITY

In Vitro

DMSO : 46 mg/mL (99.03 mM; Need ultrasonic and warming)

Concentration	Solvent	Mass		
		1 mg	5 mg	10 mg
1 mM		2.1528 mL	10.7640 mL	21.5281 mL
5 mM		0.4306 mL	2.1528 mL	4.3056 mL
10 mM		0.2153 mL	1.0764 mL	2.1528 mL

Please refer to the solubility information to select the appropriate solvent.

BIOLOGICAL ACTIVITY

Description

GKA50 is a potent glucokinase activator (EC₅₀=33 nM at 5 mM glucose). GKA50 stimulates insulin release from mouse islets of Langerhans. GKA50 is a glucose-like activator of beta-cell metabolism in rodent and human islets and a Ca²⁺-dependent modulator of insulin secretion. GKA50 shows significant glucose lowering in high fat fed female rats^{[1][2][3][4]}.

In Vitro

GKA50 (0.01-100 μM; 24 hours) enhances INS-1 cell proliferation with EC₅₀ values ranging from 1 to 2 μM^[2]. GKA50 (1.2 μM+40 μM glucose; 2-4 days) treatment reduces apoptosis induced by chronic high glucose in INS-1 cells^[2]. GKA50 activates human glucokinase enzymatic activity with an EC₅₀ of 0.022 μM. GKA50 stimulates insulin secretion in the pancreatic insulinoma cell line, INS-1, with an EC₅₀ of 0.065 μM. GKA50 reduces chronic-high-glucose-induced apoptosis via modulation of glucokinase and apoptotic protein BAD^[2].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

Cell Proliferation Assay^[2]

Cell Line: INS-1 cells (starved overnight with 3 μM glucose)

Concentration: 0.01-100 μM

	Incubation Time:	24 hours
	Result:	Stimulated cell proliferation in a dose-dependent manner, with EC ₅₀ values ranging from 1 to 2 μM.
In Vivo	GKA50 (1-30 mg/kg; p.o.) gives significant glucose lowering in an oral glucose tolerance test ^[1] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.	
	Animal Model:	High-fat-fed obese female Zucker (fa/fa) rats ^[1]
	Dosage:	1, 3, 10, 30 mg/kg
	Administration:	Oral administration
	Result:	Significant percentage glucose lowering.

REFERENCES

- [1]. Coope GJ, et al. Predictive blood glucose lowering efficacy by Glucokinase activators in high fat fed female Zucker rats. Br J Pharmacol. 2006 Oct;149(3):328-35.
- [2]. McGlasson L, et al. The glucokinase activator GKA50 causes an increase in cell volume and activation of volume-regulated anion channels in rat pancreatic β-cells. Mol Cell Endocrinol. 2011 Aug 6;342(1-2):48-53.
- [3]. Johnson D, et al. Glucose-dependent modulation of insulin secretion and intracellular calcium ions by GKA50, a glucokinase activator. Diabetes. 2007;56(6):1694-1702.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA