


MS Spectrum



Data and Time : 2017/10/15
 Sample : Enfuvirtide Ac-YF-36-NH2
 Inj. Volume : 1.000
 Mw : 4491.98
 Catalog No : HY-P0052 Batch#26748

Probe : ESI
 Nebulizer Gas Flow : 1.5L/min
 CDL : -20.0v
 CDL Temp : 250°C
 Block Temp : 400°C
 Probe bias : -4.5kv
 Detector : 1.2kv
 T.Flow : 0.2ml/min
 B.conc : 50%H2O/50%ACN