


MS Spectrum



Data and Time :2017/9/23
Sample : Amyloid beta-peptide(25-35) GM-11
Inj. Volume : 1.000
Mw : 1060.29

Probe :ESI
Nebulizer Gas Flow :1.5L/min
CDL : -20.0v
CDL Temp :250°C
Block Temp :400°C
Probe bias : -4.5kv
Detector :1.2kv
T.Flow :0.2ml/min
B.conc :50%H2O/50%ACN

Catalog No: HY-P0128 Batch#26746