

Biotinyl-(εAhx)-Gly-Arg-Gly-Asp-Ser

| | |
|----------------------|-------------------------------------------------------------------------------------------|
| Cat. No.: | HY-P4945 |
| CAS No.: | 2022956-45-0 |
| Molecular Formula: | C ₃₃ H ₅₅ N ₁₁ O ₁₂ S |
| Molecular Weight: | 829.92 |
| Sequence: | Biotinyl-{εAhx}-Gly-Arg-Gly-Asp-Ser |
| Sequence Shortening: | Biotinyl-{εAhx}-GRGDS |
| Target: | Others |
| Pathway: | Others |
| Storage: | Please store the product under the recommended conditions in the Certificate of Analysis. |

BIOLOGICAL ACTIVITY

Description

Biotinyl-(εAhx)-Gly-Arg-Gly-Asp-Ser is a polypeptide that can be found by peptide screening. Peptide screening is a research tool that pools active peptides primarily by immunoassay. Peptide screening can be used for protein interaction, functional analysis, epitope screening, especially in the field of agent research and development^[1].

REFERENCES

[1]. Birnbaum S, et al. Peptide screening. Current Opinion in Biotechnology, 1992, 3(1): 49-54.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA