

## Testican 3/SPOCK3 Protein, Mouse (HEK293, His)

Cat. No.:	HY-P71352
Synonyms:	Testican-3; Spock3
Species:	Mouse
Source:	HEK293
Accession:	Q8BKV0-1 (A23-I436)
Gene ID:	72902
Molecular Weight:	58-60 kDa

### PROPERTIES

AA Sequence	<pre> A A A V A V A G G R   S D G G N F L D E K   Q W L T T I S Q Y D   K E V G Q W N K F R D E V E D D Y F R T   W N P G K P F D Q A   L D P A K D P C L K   T K C S R H K V C I T Q D A Q T A L C I   S H R R L T H S M K   E V G G S H K Q W R   G L P S S T C K P C P I A Y A S P V C G   S D G H S Y S S Q C   K L E Y Q A C V L G   K Q I S I K C E G R C P C P S D K S M N   I G R N V K R A C S   D L E F R E V A N R   L R D W F K A L H E S G S Q N K K T K A   L L R P E R S R F D   T S I L P I C K D S   L G W M F N R L D T N Y D L L L D Q S E   L G S I Y L D K N E   Q C T K A F F N S C   D T Y K D S L I S N N E W C Y C F Q R Q   Q D P P C H T E L S   N I Q K R Q G I K K   L L G Q Y I P L C D E D G Y Y K P T Q C   H G S V G Q C W C V   D R Y G N E V V G S   R I N G V A D C A I D F E I S G D F A S   G D F R E W T D D E   G E E D D I M N D K   D D I E D D D E D E G D D D D D G D V H   D G Y I                 </pre>
Appearance	Lyophilized powder.
Formulation	Lyophilized from a 0.2 µm filtered solution of 20 mM PB, 150 mM NaCl, pH 7.4.
Endotoxin Level	<1 EU/µg, determined by LAL method.
Reconstitution	It is not recommended to reconstitute to a concentration less than 100 µg/mL in ddH <sub>2</sub> O. For long term storage it is recommended to add a carrier protein (0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose).
Storage & Stability	Stored at -20°C for 2 years. After reconstitution, it is stable at 4°C for 1 week or -20°C for longer. It is recommended to freeze aliquots at -20°C or -80°C for extended storage.
Shipping	Room temperature in continental US; may vary elsewhere.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: [tech@MedChemExpress.com](mailto:tech@MedChemExpress.com)

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA