

Screening Libraries

Proteins

Product Data Sheet

Zika virus E/Envelope Protein (107a.a, sf9, His)

Cat. No.: HY-P74459

Zika virus (ZIKV) (strain Zika SPH2015) ZIKV-E/Envelope protein (Domain III, His) Synonyms:

Species:

Sf9 insect cells Source:

Accession: ALU33341 (V593-K699)

Gene ID:

Molecular Weight: Approximately 13 kDa

PROPERTIES

Appearance	Lyophilized powder.
Formulation	Lyophilized from a 0.2 μm filtered solution of 20 mM Tris, 500 mM NaCl, pH 8.0, 10% Glycerol. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization.
Endotoxin Level	<1 EU/µg, determined by LAL method.
Reconsititution	It is not recommended to reconstitute to a concentration less than 100 $\mu g/mL$ in ddH $_2$ O.
Storage & Stability	Stored at -20°C for 2 years. After reconstitution, it is stable at 4°C for 1 week or -20°C for longer (with carrier protein). It is recommended to freeze aliquots at -20°C or -80°C for extended storage.
Shipping	Room temperature in continental US; may vary elsewhere.

DESCRIPTION

Background

Genome polyprotein is a series of protein units with similar or different functions that have been widely utilized by singlecelled or multi-cellular organisms as concentrators of countless molecular activities. Genome polyprotein is a small protein chain that is covalently linked, and it is a common means of organizing the protein set of viruses (including HIV) in nature. As the signal peptide of NS4B, genome polyprotein is essential for the anti-interferon activity of NS4B. Genome polyprotein inhibits RNA silencing by interfering with host Dicer. Genome polyprotein may play a role in viral budding. Genome polyprotein exerts cytotoxic effects by activating the mitochondrial apoptosis pathway through the M ectodomain. Genome polyprotein may display viral protein activity [1][2].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898 Fax: 609-228-5909 E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA

Page 1 of 1