

## PRELP Protein, Human (HEK293, Fc)

Cat. No.:	HY-P75981
Synonyms:	Prolargin; Proline-arginine-rich end leucine-rich repeat protein; PRELP; SLRR2A
Species:	Human
Source:	HEK293
Accession:	P51888 (Q21-I382)
Gene ID:	5549
Molecular Weight:	80-90 kDa

### PROPERTIES

Appearance	Lyophilized powder.
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization.
Endotoxin Level	<1 EU/µg, determined by LAL method.
Reconstitution	It is not recommended to reconstitute to a concentration less than 100 µg/mL in ddH <sub>2</sub> O.
Storage & Stability	Stored at -20°C for 2 years. After reconstitution, it is stable at 4°C for 1 week or -20°C for longer (with carrier protein). It is recommended to freeze aliquots at -20°C or -80°C for extended storage.
Shipping	Room temperature in continental US; may vary elsewhere.

### DESCRIPTION

Background	PRELP Protein emerges as a potential anchoring factor linking basement membranes to the underlying connective tissue, suggesting a crucial role in structural integrity. This function implies its involvement in maintaining the stability and organization of extracellular matrices. PRELP specifically interacts with the basement membrane heparan sulfate proteoglycan perlecan and triple helical collagens type I and type II, underscoring its versatile binding capabilities and its significance in mediating molecular connections within the extracellular environment. The protein's dual ability to anchor to perlecan and collagens hints at its role in orchestrating the intricate composition of basement membranes, reflecting its potential contributions to tissue architecture and homeostasis. Further exploration of PRELP's interactions and functions may unveil its specific mechanisms in basement membrane anchoring and provide insights into its broader implications in connective tissue biology.
------------	---

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: [tech@MedChemExpress.com](mailto:tech@MedChemExpress.com)

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA