

## Gag-Pol polyprotein, HIV (440a.a, His)

<b>Cat. No.:</b>	HY-P76972
<b>Synonyms:</b>	HIV-p51 / RT-p51 (group M, subtype B (isolate HXB2) Gag-Pol polyprotein Protein (His)
<b>Species:</b>	Virus
<b>Source:</b>	E. coli
<b>Accession:</b>	P04585 (P588-F1027)
<b>Gene ID:</b>	155348
<b>Molecular Weight:</b>	Approximately 54 kDa.

### PROPERTIES

<b>AA Sequence</b>	<pre> P I S P I E T V P V   K L K P G M D G P K   V K Q W P L T E E K   I K A L V E I C T E M E K E G K I S K I   G P E N P Y N T P V   F A I K K K D S T K   W R K L V D F R E L N K R T Q D F W E V   Q L G I P H P A G L   K K K K S V T V L D   V G D A Y F S V P L D E D F R K Y T A F   T I P S I N N E T P   G I R Y Q Y N V L P   Q G W K G S P A I F Q S S M T K I L E P   F R K Q N P D I V I   Y Q Y M D D L Y V G   S D L E I G Q H R T K I E E L R Q H L L   R W G L T T P D K K   H Q K E P P F L W M   G Y E L H P D K W T V Q P I V L P E K D   S W T V N D I Q K L   V G K L N W A S Q I   Y P G I K V R Q L C K L L R G T K A L T   E V I P L T E E A E   L E L A E N R E I L   K E P V H G V Y Y D P S K D L I A E I Q   K Q G Q G Q W T Y Q   I Y Q E P F K N L K   T G K Y A R M R G A H T N D V K Q L T E   A V Q K I T T E S I   V I W G K T P K F K   L P I Q K E T W E T W W T E Y W Q A T W   I P E W E F V N T P   P L V K L W Y Q L E   K E P I V G A E T F </pre>
<b>Biological Activity</b>	Data is not available.
<b>Appearance</b>	Lyophilized powder.
<b>Formulation</b>	Lyophilized from a 0.2 µm filtered solution of 50 mM Tris, 250 mM NaCl, 1 mM EDTA, 50 % glycerol, pH 7.0 or 50 mM Tris-HCL, 300 mM NaCl, pH 7.4.
<b>Endotoxin Level</b>	<1 EU/µg, determined by LAL method.
<b>Reconstitution</b>	It is not recommended to reconstitute to a concentration less than 100 µg/mL in ddH <sub>2</sub> O. For long term storage it is recommended to add a carrier protein (0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose).
<b>Storage &amp; Stability</b>	Stored at -20°C for 2 years. After reconstitution, it is stable at 4°C for 1 week or -20°C for longer (with carrier protein). It is recommended to freeze aliquots at -20°C or -80°C for extended storage.
<b>Shipping</b>	Room temperature in continental US; may vary elsewhere.

## DESCRIPTION

### Background

The Gag-Pol polyprotein assumes a central role in the virion assembly process, collaborating with the Gag polyprotein to orchestrate essential events such as binding to the plasma membrane, facilitating protein-protein interactions crucial for spherical particle formation, recruiting viral Env proteins, and packaging genomic RNA through direct interactions with the RNA packaging sequence (Psi). This polyprotein potentially regulates its own translation by binding genomic RNA in the 5'-UTR, exhibiting a dual role in promoting translation at low concentrations and encapsidating genomic RNA to inhibit translation at higher concentrations. The multipartite membrane-binding signal, including the myristoylated N-terminus, targets the polyprotein to the plasma membrane. Additionally, the Matrix protein within the polyprotein is implicated in the pre-integration complex and is associated with the release from the host cell mediated by Vpu. The ability of Gag-Pol to bind to RNA further underscores its multifaceted functions in the intricate process of virion assembly.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: [tech@MedChemExpress.com](mailto:tech@MedChemExpress.com)

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA