## Product Data Sheet

## Inhibitors • Screening Libraries • Proteins



## hsa-miR-5004-3p mimic

| Cat. No.:         | HY-R01493   |
|-------------------|---|
| Molecular Weight: | 14035.5   |
| Target:           | MicroRNA  |
| Pathway:          | Epigenetics   |
| Storage:          | Please store the product under the recommended conditions in the Certificate of Analysis. |

| BIOLOGICAL ACTIVITY |  |  |
|---------------------|--|--|
| Description         | hsa-miR-5004-3p mimics are small, chemically synthesized double-stranded RNAs that mimic endogenous miRNAs and<br>enable miRNA functional analysis by up-regulation of miRNA activity.   |  |
| In Vitro            | <ol> <li>miRNA Resuspension</li> <li>Briefly centrifuge the tube to ensure that the dried miRNA is at the bottom of the tube.</li> <li>Resuspend the miRNA using nuclease free water to generate 20 µM stock solution.</li> <li>For 5 nmol miRNA: add 250 µL nuclease free water.</li> <li>For 20 nmol miRNA: add 1000 µL nuclease free water.</li> <li>A liquot miRNA: add 1000 µL nuclease free water.</li> <li>A liquot miRNA: add 20° µ nuclease free water.</li> <li>A liquot miRNA: add 20° µ nuclease free water.</li> <li>A liquot miRNA: add 000 µL nuclease free water.</li> <li>A liquot miRNA: add 000 µL nuclease free water.</li> <li>A liquot miRNA: add 1000 µL nuclease free water.</li> <li>A liquot miRNA: add 1000 µL nuclease free water.</li> <li>A liquot miRNA: add 20° µ nuclease free water.</li> <li>A liquot miRNA: add 1000 µL nuclease free water.</li> <li>A liquot miRNA: add 1000 µL nuclease free water.</li> <li>A liquot miRNA: add 1000 µL nuclease free water.</li> <li>A liquot miRNA: add 1000 µL nuclease free water.</li> <li>A liquot miRNA: add 20° µL nuclease free water.</li> <li>A liquot miRNA: add 1000 µL nuclease free water.</li> <li>A liquot miRNA: add 1000 µL nuclease free water.</li> <li>A liquot miRNA: add 000 µL nuclease free water.</li> <li>A liquot miRNA: add 200 µL nuclease free water.</li> <li>A liquot miRNA: add 200 µL nuclease free water.</li> <li>A liquot miRNA: add 200 µL nuclease free water.</li> <li>A fransfection result.</li> <li>Transfection mix A and B.</li> <li>For per well of a 6-well plate: A: 120 µL serum-free medium + 5 µL miRNA mimics; B: 121 µL serum-free medium + 4 µL PolyFast Transfection Reagent (HY-K1014).</li> <li>For per well of a 24-well plate: A: 23.75 µL serum-free medium + 1.25 µL miRNA mimics; B: 24 µL serum-free medium + 1 µL PolyFast Transfection Reagent (HY-K1014).</li> <li>For per well of a 96-well plate: A: 4.75 µL serum-free medium + 0.25 µL miRNA mimics; B: 4.</li></ol> |  |

optimization experiments using varying mimic/inhibitor concentrations should be performed. The optimized range suggests changing the miRNA concentration in the range of 10 to 200nM.

If other transfection reagents are used, the amount of transfection reagent needs to be adjusted according to the specific situation.

3.2 Mix A and B gently. Incubate at room temperature for 15 minutes.

3.3 Remove culture medium from cells, wash with PBS.

3.4 Add transfection mix (A+B) to cells.

For per well of a 6-well plate: add 1750  $\mu$ L fresh medium without Pen/Strep, then add 250  $\mu$ L of the transfection mix (A+B) to the well, and mix well.

For per well of a 24-well plate: add 450  $\mu$ L fresh medium without Pen/Strep, then add 50  $\mu$ L of the transfection mix (A+B) to the well, and mix well.

For per well of a 96-well plate: add 90  $\mu$ L fresh medium without Pen/Strep, then add 10  $\mu$ L of the transfection mix (A+B) to the well, and mix well.

3.5 Incubate cells for 1–3 days at 37°C. Then, analyze transfected cells. The medium can be replaced with fresh serumcontaining medium after 6 hours if necessary.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898 Fax: 609-228-5909 E-mail: tech@MedChemExpress.com Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA