

## MDM2 Human Pre-designed siRNA Set A

|                  |   |
|------------------|---|
| <b>Cat. No.:</b> | HY-RS08264  |
| <b>Target:</b>   | Small Interfering RNA (siRNA)   |
| <b>Pathway:</b>  | Epigenetics   |
| <b>Storage:</b>  | Please store the product under the recommended conditions in the Certificate of Analysis. |

### BIOLOGICAL ACTIVITY

#### Description

MDM2 Human Pre-designed siRNA Set A contains three designed siRNAs for MDM2 gene (Human), as well as a negative control, a positive control, and a FAM-labeled negative control.

#### Components

MDM2 siRNA-1: 5 nmol (HPLC)

MDM2 siRNA-2: 5 nmol (HPLC)

MDM2 siRNA-3: 5 nmol (HPLC)

siRNA Negative Control: 5 nmol (HPLC)

FAM-labeled siRNA Negative Control: 5 nmol (HPLC)

GAPDH siRNA Positive Control: 5 nmol (HPLC)

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: [tech@MedChemExpress.com](mailto:tech@MedChemExpress.com)

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA