Inhibitors

tans-4-Hydroxy-D-proline hydrochloride

Cat. No.: HY-W003511

CAS No.: 142347-81-7

Molecular Formula: $C_5H_{10}CINO_3$ Molecular Weight: 167.59

Target: ADC Linker; PROTAC Linkers

Pathway: Antibody-drug Conjugate/ADC Related; PROTAC

Storage: Powder -20°C 3 years

 $\begin{array}{ccc} & 4^{\circ}\text{C} & 2 \text{ years} \\ \text{In solvent} & -80^{\circ}\text{C} & 6 \text{ months} \\ & -20^{\circ}\text{C} & 1 \text{ month} \end{array}$

HCI

BIOLOGICAL ACTIVITY

Description	tans-4-Hydroxy-D-proline hydrochloride is a non-cleavable ADC linker used in the synthesis of antibody-drug conjugates (ADCs). tans-4-Hydroxy-D-proline hydrochloride is also a alkyl chain-based PROTAC linker that can be used in the synthesis of PR
IC ₅₀ & Target	Non-cleavable Linker
In Vitro	ADCs are comprised of an antibody to which is attached an ADC cytotoxin through an ADC linker ^[1] . PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins ^[2] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.

REFERENCES

[1]. Beck A, et al. Strategies and challenges for the next generation of antibody-drug conjugates. Nat Rev Drug Discov. 2017;16(5):315-337.

[2]. Nalawansha DA, et al. PROTACs: An Emerging Therapeutic Modality in Precision Medicine. Cell Chem Biol. 2020;27(8):998-985.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA