Fmoc-Thr(tBu)-OH

MedChemExpress

Cat. No.:	HY-W00770	6		
CAS No.:	71989-35-0			
Molecular Formula:	C ₂₃ H ₂₇ NO ₅			
Molecular Weight:	397.46			
Target:	Amino Acid Derivatives			
Pathway:	Others			
Storage:	Powder	-20°C	3 years	
		4°C	2 years	
	In solvent	-80°C	6 months	
		-20°C	1 month	

®

SOLVENT & SOLUBILITY

	Preparing Stock Solutions	Solvent Mass Concentration	1 mg	5 mg	10 mg		
		1 mM	2.5160 mL	12.5799 mL	25.1598 mL		
		5 mM	0.5032 mL	2.5160 mL	5.0320 mL		
		10 mM	0.2516 mL	1.2580 mL	2.5160 mL		
	Please refer to the so	lubility information to select the app	propriate solvent.				
In Vivo	1. Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline Solubility: 2.5 mg/mL (6.29 mM); Suspended solution; Need ultrasonic						
	2. Add each solvent one by one: 10% DMSO >> 90% (20% SBE-β-CD in saline) Solubility: ≥ 2.5 mg/mL (6.29 mM); Clear solution						
	3. Add each solvent one by one: 10% DMSO >> 90% corn oil Solubility: ≥ 2.5 mg/mL (6.29 mM); Clear solution						

BIOLOGICAL ACTIVITY				
Description	Fmoc-Thr(tBu)-OH is a threonine derivative ^[1] .			
In Vitro	Amino acids and amino acid derivatives have been commercially used as ergogenic supplements. They influence the secretion of anabolic hormones, supply of fuel during exercise, mental performance during stress related tasks and prevent exercise induced muscle damage. They are recognized to be beneficial as ergogenic dietary substances ^[1] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.			

Y^HN O 0[∕]

ОН

REFERENCES

[1]. Luckose F, et al. Effects of amino acid derivatives on physical, mental, and physiological activities. Crit Rev Food Sci Nutr. 2015;55(13):1793-1144.

Caution: Product has not been fully validated for medical applications. For research use only.

 Tel: 609-228-6898
 Fax: 609-228-5909
 E-mail: tech@MedChemExpress.com

 Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA