

# **Product** Data Sheet


# Polyvinyl alcohol (Mw 146000-186000, 87-89% hydrolyzed)

Cat. No.:HY-Y0850PCAS No.:9002-89-5Molecular Formula: $(C_2H_4O)_X$ Target:Others

Storage: Please store the product under the recommended conditions in the Certificate of

Analysis.

Others



## **BIOLOGICAL ACTIVITY**

### Description

Pathway:

Polyvinyl alcohol (Mw 146000-186000, 87-89% hydrolyzed) is a polyvinyl alcohol with a molecular weight of 146000-186000 and hydrolytic properties. The degree of hydrolysis refers to the degree to which the acetate groups in the original polyvinyl acetate are converted into hydroxyl groups during the hydrolysis process. Polyvinyl alcohol (Mw 146000-186000, 87-89% hydrolyzed) is the hydrolysis and removal of acetate groups after the polymerization of ethylene acetate. And polyvinyl alcohol is obtained. A degree of hydrolysis of 87-89% indicates that a large part of the acetate groups have been removed, resulting in a large number of hydroxyl groups in the PVA structure. Polyvinyl alcohol with different degrees of hydrolysis can be used to self-crosslink to form cryogel, which can be used as biological excipient<sup>[1]</sup>.

### **REFERENCES**

[1]. Elsherbiny DA, et al. Self-crosslinked polyvinyl alcohol/cellulose nanofibril cryogels loaded with synthesized aminophosphonates as antimicrobial wound dressings. J Mater Chem B. 2023 Aug 2;11(30):7144-7159.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA

Inhibitors • Sc

Screening Libraries

**Proteins**