

Affinity Chromatography Column

1 Contents


Cat. No.	Product Name	Package
HY-K0221-1 mL-10 pcs	Affinity Chromatography Column-1 mL	10 pcs
HY-K0221-3 mL-10 pcs	Affinity Chromatography Column-3 mL	10 pcs
HY-K0221-6 mL-10 pcs	Affinity Chromatography Column-6 mL	10 pcs
HY-K0221-12 mL-10 pcs	Affinity Chromatography Column-12 mL	10 pcs
HY-K0221-30 mL-5 pcs	Affinity Chromatography Column-30 mL	5 pcs
HY-K0221-60 mL-5 pcs	Affinity Chromatography Column-60 mL	5 pcs

2 General Information

MCE Affinity Chromatography (AC) Column offers empty columns perfectly designed for laboratory applications of purification for recombinant proteins with different tags, enzymes, antibodies and antigens.

MCE Affinity Chromatography Column is made of high-purity polypropylene with a very low adsorption to various biomolecules. Frits for AC column are made of UHMWPE, which is hydrophilic, and optimized for minimizing non-specific adsorption.

3 Column Reference




4 Application

Load of AC media: Transfer an appropriate-sized affinity chromatography media to the empty AC column with packed between two frits.

The operation after loading can be seen in the relevant media product instructions.

Procedure Reference:



MedChemExpress USA

Tel: 609-228-6898 Fax: 609-228-5909

E-mail: sales@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA